

Product datasheet

Anti-DDX41 antibody [2254D3] ab50704

画像数 1

製品の概要

製品名	Anti-DDX41 antibody [2254D3]
製品の詳細	Mouse monoclonal [2254D3] to DDX41
由来種	Mouse
アプリケーション	適用あり: WB, Dot blot
種交差性	交差種: Human
免疫原	Recombinant human DDX41 fragment

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	Preservative: 0.05% Sodium Azide Constituents: 1% BSA, PBS (8mM Sodium phosphate, 3mM Potassium chloride, 140mM Sodium chloride, 1.5mM Potassium phosphate, pH 7.4).
精製度	Protein G purified
特記事項(精製)	Ab50704 was purified using protein G column chromatography from culture supernatant of hybridoma cultured in a medium containing bovine IgG-depleted (approximately 95%) fetal bovine serum.
ポリ/モノ	モノクローナル
クローン名	2254D3
アイソタイプ	IgG1

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab50704** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent dilution. Predicted molecular weight: 69 kDa.

アプリケーション	Abreviews	特記事項
Dot blot		Use at an assay dependent dilution.

## ターゲット情報

機能	Probable ATP-dependent RNA helicase. Is required during post-transcriptional gene expression. May be involved in pre-mRNA splicing.
配列類似性	Belongs to the DEAD box helicase family, DDX41 subfamily. Contains 1 CCHC-type zinc finger. Contains 1 helicase ATP-binding domain. Contains 1 helicase C-terminal domain.
細胞内局在	Nucleus.

## 画像



Anti-DDX41 antibody [2254D3] (ab50704) +  
Recombinant DDX41

**Predicted band size:** 69 kDa

**Observed band size:** 40 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors