



Product datasheet

Anti-DMRTB1 antibody ab30911

画像数 1

製品の概要

製品名	Anti-DMRTB1 antibody
製品の詳細	Rabbit polyclonal to DMRTB1
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Synthetic peptide within Human DMRTB1 aa 222-271 (internal sequence). The exact sequence is proprietary. Sequence: LGYLDAPPGVPLQQGFRHVSRSQYQGGGLVSEPGGDFQPSYYLPPPPPL Database link: Q96MA1
	 Run BLAST with  Run BLAST with

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: None Constituents: 2% Sucrose, PBS.
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab30911** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリ ケーショ ン	Abreviews	特記事項
------------------	-----------	------

WB Use a concentration of 2.5 µg/ml. Detects a band of approximately 44 kDa (predicted molecular weight: 36 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

ターゲット情報

組織特異性	Testis.
配列類似性	Belongs to the DMRT family. Contains 1 DM DNA-binding domain.
細胞内局在	Nucleus.

画像



Anti-DMRTB1 antibody (ab30911) at 2.5 µg/ml + HepG2 Cell Lysate at 10 µg

Secondary

HRP conjugated anti rabbit IgG at 1/50,000 - 1/100,000

Predicted band size: 36 kDa

Observed band size: 44 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors