

Product datasheet

Anti-E tag antibody [10B11] ab106184

画像数 1

製品の概要

<b>製品名</b>	Anti-E tag antibody [10B11]
<b>製品の詳細</b>	Mouse monoclonal [10B11] to E tag
<b>アプリケーション</b>	<b>適用あり:</b> WB, ELISA
<b>免疫原</b>	Synthetic peptide corresponding to a fragment ( GAPVPYPDPLEPR ) E Tag, expressed in E. coli.

 [Run BLAST with](#)

 [Run BLAST with](#)

**ポジティブ・コントロール** E tagged fusion protein

**特記事項** Abcam is committed to meeting high standards of ethical manufacturing and as such, we will be discontinuing this product, which has been generated by the ascites method, within the next year. We are sorry for any inconvenience this may cause. We would recommend antibody [ab3397](#) as a replacement.

製品の特性

<b>製品の状態</b>	Liquid
<b>保存方法</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
<b>バッファー</b>	Preservative: None Constituents: dH2O
<b>精製度</b>	Protein G purified
<b>ポリ/モノ</b>	モノクローナル
<b>クローン名</b>	10B11
<b>アイソタイプ</b>	IgG1

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab106184** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB		1/1000.
ELISA		Use at an assay dependent dilution.

## ターゲット情報

### 関連性

E tag is a peptide epitope tag, sequence GAPVPYPDPLEPR.

## 画像



Anti-E tag antibody [10B11] (ab106184) at  
1/1000 dilution + E tagged fusion protein

Western blot - E tag antibody [10B11] (ab106184)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors