

Product datasheet

Anti-E tag antibody [10B11] ab106184

画像数 1

製品の概要

製品名	Anti-E tag antibody [10B11]
製品の詳細	Mouse monoclonal [10B11] to E tag
由来種	Mouse
アプリケーション	適用あり: WB, ELISA
免疫原	Synthetic peptide corresponding to a fragment (GAPVPYPDPLEPR) E Tag, expressed in E. coli. Run BLAST with Run BLAST with
ポジティブ・コントロール	E tagged fusion protein
特記事項	Abcam is committed to meeting high standards of ethical manufacturing and as such, we will be discontinuing this product, which has been generated by the ascites method, within the next year. We are sorry for any inconvenience this may cause. We would recommend antibody ab3397 as a replacement.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
バッファー	Preservative: None Constituents: dH2O
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	10B11
アイソタイプ	IgG1

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab106184** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB		1/1000.
ELISA		Use at an assay dependent dilution.

ターゲット情報

関連性

E tag is a peptide epitope tag, sequence GAPVPYPDPLEPR.

画像



Anti-E tag antibody [10B11] (ab106184) at
1/1000 dilution + E tagged fusion protein

Western blot - E tag antibody [10B11] (ab106184)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors