

Product datasheet

Anti-ELK4 antibody ab49244

画像数 1

製品の概要

製品名	Anti-ELK4 antibody
製品の詳細	Rabbit polyclonal to ELK4
由来種	Rabbit
アプリケーション	適用あり: WB, ELISA
種交差性	交差種: Human
免疫原	A region within synthetic peptide: VAPLSPARLQ GANTLFQFPS VLNSHGPFTL SGLDGPSTPG PFSPDLQKT, corresponding to C terminal amino acids 383-431 of Human ELK4 Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	Jurkat cell lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: None Constituents: 2% Sucrose, PBS
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab49244** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 2.5 µg/ml. Detects a band of approximately 41 kDa (predicted molecular weight: 47 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.


アプリ ケーショ ン	Abreviews	特記事項
------------------	-----------	------

ELISA Use at an assay dependent concentration.
ELISA titre using peptide based assay: 1/62500.

ターゲット情報

機能	Forms a ternary complex with the serum response factor (SRF). Requires DNA-bound SRF for ternary complex formation and makes extensive DNA contacts to the 5' side of SRF, but does not bind DNA autonomously.
配列類似性	Belongs to the ETS family. Contains 1 ETS DNA-binding domain.
細胞内局在	Nucleus.

画像



Anti-ELK4 antibody (ab49244) at 2.5 µg/ml + Jurkat cell lysate at 10 µg with 5% skim milk in PBS.

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000 dilution

Predicted band size: 47 kDa

Observed band size: 41 kDa

Gel concentration: 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors