

Product datasheet

Anti-EME1 antibody ab88878

★☆☆☆☆ 2 Abreviews 画像数 3

製品の概要

製品名	Anti-EME1 antibody
製品の詳細	Mouse polyclonal to EME1
由来種	Mouse
アプリケーション	適用あり: WB, ICC/IF
種交差性	交差種: Human
免疫原	Full length protein, corresponding to amino acids 1 - 583 of Human EME1 (AAH16470.1)
ポジティブ・コントロール	HeLa or Jurkat cell lysate. HeLa cells.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: 1X PBS, pH 7.2
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab88878** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB	★☆☆☆☆	Use a concentration of 1 µg/ml. Predicted molecular weight: 63 kDa.
ICC/IF		Use a concentration of 10 µg/ml.

ターゲット情報


機能	Interacts with MUS81 to form a DNA structure-specific endonuclease with substrate preference for branched DNA structures with a 5'-end at the branch nick. Typical substrates include 3'-flap structures, replication forks and nicked Holliday junctions. May be required in mitosis for the processing of stalled or collapsed replication forks.
配列類似性	Belongs to the EME1/MMS4 family.
細胞内局在	Nucleus > nucleolus. Recruited to regions of DNA damage in S-phase cells.

画像



Anti-EME1 antibody (ab88878) at 1 µg/ml +
HeLa cell lysate at 50 µg

Predicted band size: 63 kDa



Anti-EME1 antibody (ab88878) at 1 µg/ml +
Jurkat cell lysate
at 50 µg

Predicted band size: 63 kDa



ab88878 at 10 µg/ml staining EME1 in HeLa
cells by immunofluorescence.

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors