

## Product datasheet

### EN2 peptide (320-333) ab45879

#### 製品の概要

---

製品名 EN2 peptide (320-333)

#### 製品の詳細

---

由来 Synthetic

#### アミノ酸配列

生物種 Human

配列 C-NHSTTAKEGKSDSE

領域 320 to 333

#### 特性

---

Our [Abpromise guarantee](#) covers the use of **ab45879** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション Blocking - Blocking peptide for Anti-EN2 antibody ([ab45867](#))

製品の状態 Liquid

#### 前処理および保存

---

保存方法および安定性 Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.

#### 関連情報

---

関連性 EN1 and EN2 encode homeodomain containing proteins and have been implicated in the control of pattern formation during development of the central nervous system. EN2 may play a role in susceptibility to autism spectrum disorders.

細胞内局在 Nuclear

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors