

## Product datasheet

# EN2 peptide (320-333) ab45879

### 製品の概要

**製品名** EN2 peptide (320-333)

### 製品の詳細

**由来** Synthetic

#### アミノ酸配列

**生物種** Human

**配列** C-NHSTTAKEGKSDSE

**領域** 320 to 333

### 特性

Our [Abpromise guarantee](#) covers the use of **ab45879** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**アプリケーション** Blocking - Blocking peptide for Anti-EN2 antibody ([ab45867](#))

**製品の状態** Liquid

### 前処理および保存

**保存方法および安定性** Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.

### 関連情報

**関連性** EN1 and EN2 encode homeodomain containing proteins and have been implicated in the control of pattern formation during development of the central nervous system. EN2 may play a role in susceptibility to autism spectrum disorders.

**細胞内局在** Nuclear

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors