

Product datasheet

Anti-EVI5 antibody ab70790

★★★★★ 2 Abreviews [画像数 1](#)

製品の概要

| | |
|--------------|--|
| 製品名 | Anti-EVI5 antibody |
| 製品の詳細 | Rabbit polyclonal to EVI5 |
| 由来種 | Rabbit |
| アプリケーション | 適用あり: IP 適用なし: WB |
| 種交差性 | 交差種: Human 交差が予測される動物種: Horse, Guinea pig, Cow, Dog, Pig, Chimpanzee, Rhesus monkey, Gorilla, Orangutan, Elephant |
| 免疫原 | A region between residue 525 and 575 of human EVI5. (NP_005656.4) |
| ポジティブ・コントロール | HeLa whole cell lysate |

製品の特性

| | |
|--------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |
| バッファー | Preservative: 0.09% Sodium Azide Constituents: 8mM PBS, 60mM Citrate, 150mM Tris, pH 7-8 |
| 精製度 | Immunogen affinity purified |
| ポリ/モノ | ポリクローナル |
| アイソタイプ | IgG |

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab70790** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|-----------------------------|
| IP | | Use at 2-5 µg/mg of lysate. |

追加情報 Is unsuitable for WB.

ターゲット情報

| | |
|-------|--|
| 関連性 | EVI5 is a novel centrosomal protein with a complex expression pattern and subcellular localization, possibly involved in centrosome stability and dynamics. It may act as an oncogene. |
| 細胞内局在 | Nucleus, centrosome. |

画像



HeLa whole cell lysate (1 mg for IP, 20% of IP loaded) with ab70790 at 3 μ g/mg lysate for IP.
Detection: Chemiluminescence with an exposure time of 1 second.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors