


## Product datasheet

### Anti-FMO5 antibody ab103973

画像数 1

#### 製品の概要

製品名	Anti-FMO5 antibody
製品の詳細	Rabbit polyclonal to FMO5
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Rabbit, Cow 
免疫原	Recombinant protein fragment containing a sequence corresponding to a region within amino acids 1 and 172 of Human FMO5 (NM_001461)
ポジティブ・コントロール	HCT116 whole cell lysate; A549, HeLa and HepG2 cell lysates

#### 法規制情報

**医薬用外毒物**

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 20% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7.0
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab103973** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/3000. Predicted molecular weight: 60 kDa.

## ターゲット情報

機能	In contrast with other forms of FMO it does not seem to be a drug-metabolizing enzyme.
組織特異性	Expressed in fetal and adult liver.
配列類似性	Belongs to the FMO family.
細胞内局在	Microsome membrane. Endoplasmic reticulum membrane.

## 画像



Anti-FMO5 antibody (ab103973) at 1/1000  
dilution + HCT116 whole cell lysate  
at 30 µg

**Predicted band size:** 60 kDa

12% SDS PAGE

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors