

Product datasheet

Anti-FOLR2 antibody ab56067

[1 Abreviews](#) [2 References](#)

製品の概要

製品名	Anti-FOLR2 antibody
製品の詳細	Rabbit polyclonal to FOLR2
由来種	Rabbit
アプリケーション	適用あり: Flow Cyt, WB, ELISA
種交差性	交差種: Human, Recombinant fragment 交差が予測される動物種: Mouse
免疫原	Synthetic peptide corresponding to Human FOLR2 (N terminal). Sequence: CMDAKHHKTKPGPEDKLHD <div style="display: flex; justify-content: flex-end; gap: 20px;"> Run BLAST with Run BLAST with </div>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.01% Sodium Azide Constituents: 50% Glycerol, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab56067** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

Flow Cyt Use at an assay dependent concentration. PubMed: 23439585
[ab171870](#) - Rabbit polyclonal IgG, is suitable for use as an isotype control with this antibody.

WB Use at an assay dependent concentration. Predicted molecular weight: 29 kDa. This antibody has been tested in Western blot against the recombinant peptide used as an immunogen. We have no data on detection of endogenous protein.

ELISA Use at an assay dependent concentration.

ターゲット情報

機能	Binds to folate and reduced folic acid derivatives and mediates delivery of 5-methyltetrahydrofolate to the interior of cells.
組織特異性	Expressed in placenta and hematopoietic cells. Expression is increased in malignant tissues.
配列類似性	Belongs to the folate receptor family.
翻訳後修飾	Eight disulfide bonds are present. N-glycosylated.
細胞内局在	Cell membrane. Secreted.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors