

Product datasheet

Anti-FOXE1 antibody ab105815

画像数 1

製品の概要

製品名	Anti-FOXE1 antibody
製品の詳細	Rabbit polyclonal to FOXE1
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Mouse 交差が予測される動物種: Rat, Cow
免疫原	Synthetic peptide from within C terminal amino acids 324-372 (FYGRTSPGQF GAALGPCYNP GGQLGAGGGG AYHSRHATAY PGAVDRFVSA) of Mouse FOXE1 (NP_899121) Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	Mouse liver lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab105815** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 38 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

ターゲット情報

機能	Probable transcription factor. Could be involved in thyroid gland organogenesis.
組織特異性	Detected in adult brain, placenta, lung, liver, skeletal muscle, kidney, pancreas, heart, colon, small intestine testis and thymus. Expression was strongest in heart and pancreas.
関連疾患	Defects in FOXE1 are the cause of Bamforth-Lazarus syndrome (BLS) [MIM:241850]. BLS is associated with thyroid agenesis, cleft palate and choanal atresia.
配列類似性	Contains 1 fork-head DNA-binding domain.
翻訳後修飾	Phosphorylated.
細胞内局在	Nucleus.

画像



Anti-FOXE1 antibody (ab105815) at 1 μ g/ml +
mouse liver lysate at 10 μ g

Predicted band size: 38 kDa

Gel concentration 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors