

## Product datasheet

# Anti-Fbp5B antibody ab96858

画像数 2

### 製品の概要

製品名	Anti-Fbp5B antibody
製品の詳細	Rabbit polyclonal to Fbp5B
由来種	Rabbit
アプリケーション	<b>適用あり:</b> WB, IHC-P
種交差性	<b>交差種:</b> Human
免疫原	Synthetic peptide containing a sequence corresponding to a region within amino acids 646-708 (PAKYQPYKKR GLCSRTACGF DFCVLCLCAY HGSEEC SRGA AKPRNRKDAL PGSAQSKRNL KRL) of Human Fbp5B (NP_001025031). <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
ポジティブ・コントロール	293T whole cell lysate; Saos-2; A431, H1299, HeLa, HepG2, Molt-4 and Raji cell lysates

### 法規制情報

**医薬用外毒物**

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 10% Glycerol, 1X PBS, pH 7.0
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab96858** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/3000. Predicted molecular weight: 78 kDa.
IHC-P		1/100 - 1/250.

## ターゲット情報

**関連性** Members of the F-box protein family, such as Fbp5B, are characterized by an approximately 40-amino acid F-box motif. SCF complexes, formed by SKP1 (MIM 601434), cullin (see CUL1; MIM 603134), and F-box proteins, act as protein-ubiquitin ligases. F-box proteins interact with SKP1 through the F box, and they interact with ubiquitination targets through other protein interaction domains (Jin et al., 2004 [PubMed 15520277]).


## 画像



Anti-Fbp5B antibody (ab96858) at 1/1000 dilution + 293T whole cell lysate at 30 µg

**Predicted band size:** 78 kDa

7.5% SDS-PAGE



Immunohistochemical analysis of Fbp5B in paraffin-embedded Saos-2, using ab96858 at 1/100 dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors