

Product datasheet

Anti-GALR3 antibody ab48008

画像数 2

製品の概要

製品名	Anti-GALR3 antibody
製品の詳細	Goat polyclonal to GALR3
由来種	Goat
アプリケーション	<b>適用あり:</b> IHC-P, ELISA, WB
種交差性	<b>交差種:</b> Human
免疫原	Synthetic peptide: C-QGPEPREGPVHGGE , corresponding to C terminal amino acids 349-362 of Human GALR3. <a href="#">Run BLAST with</a> <a href="#">Run BLAST with</a>
ポジティブ・コントロール	Human skin and testes lysates.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.02% Sodium Azide Constituents: 0.5% BSA, Tris buffered saline, pH 7.3
精製度	Immunogen affinity purified
特記事項(精製)	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab48008** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P		Use a concentration of 4 - 6 µg/ml. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
ELISA		1/32000.
WB		Use a concentration of 1 - 3 µg/ml. Detects a band of approximately 37 kDa (predicted molecular weight: 40 kDa).

## ターゲット情報

機能	Receptor for the hormone galanin.
配列類似性	Belongs to the G-protein coupled receptor 1 family.
細胞内局在	Cell membrane.

## 画像




Anti-GALR3 antibody (ab48008) at 2 µg/ml +  
Human skin lysate in RIPA buffer at 35 µg/ml

Developed using the ECL technique.

**Predicted band size:** 40 kDa

**Observed band size:** 37 kDa



ab48008 (3.8µg/ml) staining of paraffin embedded Human Testis shows cellular membrane staining of spermatocytes.  
Steamed antigen retrieval with citrate buffer pH 6, AP-staining.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors