

Product datasheet

Anti-GALR3 antibody ab48008

画像数 2

製品の概要

製品名	Anti-GALR3 antibody
製品の詳細	Goat polyclonal to GALR3
アプリケーション	<b>適用あり:</b> IHC-P, ELISA, WB
種交差性	<b>交差種:</b> Human
免疫原	Synthetic peptide: C-QGPEPREGPVHGGE , corresponding to C terminal amino acids 349-362 of Human GALR3. <a href="#">Run BLAST with</a> <a href="#">Run BLAST with</a>
ポジティブ・コントロール	Human skin and testes lysates.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.02% Sodium Azide Constituents: 0.5% BSA, Tris buffered saline, pH 7.3
精製度	Immunogen affinity purified
特記事項 (精製)	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
ポリモノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab48008** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P		Use a concentration of 4 - 6 µg/ml. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
ELISA		1/32000.
WB		Use a concentration of 1 - 3 µg/ml. Detects a band of approximately 37 kDa (predicted molecular weight: 40 kDa).

## ターゲット情報

機能	Receptor for the hormone galanin.
配列類似性	Belongs to the G-protein coupled receptor 1 family.
細胞内局在	Cell membrane.


## 画像



Anti-GALR3 antibody (ab48008) at 2 µg/ml +  
Human skin lysate in RIPA buffer at 35 µg/ml  
Developed using the ECL technique

**Predicted band size** : 40 kDa

**Observed band size** : 37 kDa



ab48008 (3.8µg/ml) staining of paraffin embedded Human Testis shows cellular membrane staining of spermatocytes.  
Steamed antigen retrieval with citrate buffer pH 6, AP-staining.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors