

Product datasheet

Anti-GM-CSF antibody (Biotin) ab84259

製品の概要

製品名	Anti-GM-CSF antibody (Biotin)
製品の詳細	Rabbit polyclonal to GM-CSF (Biotin)
由来種	Rabbit
標識	Biotin
アプリケーション	<b>適用あり:</b> Sandwich ELISA, WB
種交差性	<b>交差種:</b> Human
免疫原	Highly pure (>98%) recombinant human GM-CSF.
ポジティブ・コントロール	

[ウェスタン・ブロット用ポジティブコントロール  
Recombinant human GM-CSF protein](#) >

製品の特性

製品の状態	Lyophilised:Centrifuge vial prior to opening. Reconstitute in sterile PBS containing 0.1% BSA to a concentration of 0.1-1.0 mg/ml.
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: None Constituents: PBS, pH 7.2
精製度	Protein A purified
特記事項(精製)	ab84259 is sterile filtered.
ポリモノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab84259** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Sandwich ELISA		Use a concentration of 0.25 - 1 µg/ml. Allows the detection of at least 0.2 - 0.4 ng/well of recombinant human GM-CSF. Can be paired for Sandwich ELISA with <a href="#">ab9667</a> . Can be used as detection antibody when paired with recommended pair.
WB		Use a concentration of 0.1 - 0.2 µg/ml. Predicted molecular weight: 17 kDa. The detection limit for recombinant human GM-CSF is 1.5 - 3.0 ng/lane under either reducing or non-reducing conditions.

## ターゲット情報

機能	Cytokine that stimulates the growth and differentiation of hematopoietic precursor cells from various lineages, including granulocytes, macrophages, eosinophils and erythrocytes.
配列類似性	Belongs to the GM-CSF family.
細胞内局在	Secreted.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors