

Product datasheet

Anti-GPBP1L1 antibody ab122002

画像数 2

製品の概要

製品名	Anti-GPBP1L1 antibody
製品の詳細	Rabbit polyclonal to GPBP1L1
アプリケーション	適用あり: ICC/IF, IHC-P
種交差性	交差種: Human
免疫原	Recombinant Protein Epitope Signature Tag: GMSQRSGGGT GNHRHWNGSF HSRKGCAFQE KPPMEIREEK KEDKVEKLQF EEEDFPSLNP EAGKQHQPGR PIGTPSGVWE NPPSAKQP, corresponding to amino acids 105-192 of Human GPBP1L1. Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	Human cerebellum tissue Transfected HEK293T over-expression lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 59% PBS, 40% Glycerol
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab122002** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
ICC/IF		Use a concentration of 1 - 4 µg/ml. Recommend PFA Fixation and Triton X-100 treatment

IHC-P 1/200 - 1/500. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

ターゲット情報

機能	Possible transcription factor.
配列類似性	Belongs to the vasculin family.
細胞内局在	Nucleus.


Anti-GPBP1L1 antibody 画像



Immunofluorescent staining of Human cell line A-431 shows positivity in nucleus but not nucleoli and plasma membrane.

Recommended concentration of ab122002 1-4 µg/ml. Cells treated with PFA/Triton X-100.

Immunocytochemistry/ Immunofluorescence - Anti-GPBP1L1 antibody (ab122002)



ab122002, at a 1/350 dilution, staining GPBP1L1 in paraffin-embedded Human cerebellum tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-GPBP1L1 antibody (ab122002)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors