

## Product datasheet

# Anti-GST antibody (FITC) ab106228

### 画像数 1

#### 製品の概要

製品名	Anti-GST antibody (FITC)
製品の詳細	Rabbit polyclonal to GST (FITC)
由来種	Rabbit
標識	FITC. Ex: 493nm, Em: 528nm
アプリケーション	<b>適用あり:</b> WB
免疫原	Full length GST isolated from Schistosoma japonicum.
ポジティブ・コントロール	Recombinant GST protein.

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.20 Preservative: 0.01% Sodium azide Constituents: 0.42% Potassium phosphate, 0.88% Sodium chloride, 1% BSA
精製度	Immunogen affinity purified
特記事項 (精製)	ab106228 was prepared from monospecific antiserum by immunoaffinity chromatography using GST coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities and extensive dialysis against the buffer.
ポリモノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab106228** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB		1/1000 - 1/10000. Predicted molecular weight: 26 kDa.

## ターゲット情報

### 関連性

GST (Glutathione S-Transferase) is a 26kDa protein encoded by the parasitic helminth *Schistosoma japonicum* and widely used in the pGEX family of GST plasmid expression vectors as a fusion protein with foreign proteins.

### 画像



Anti-GST antibody (FITC) (ab106228) at  
1/1000 dilution + Recombinant GST protein at  
0.25 µg

#### Secondary

HRP conjugated Goat Anti-Rabbit at 1/40000  
dilution

**Predicted band size:** 26 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors