

Product datasheet

Anti-Galectin 3 antibody ab2473

製品の概要

| | |
|----------|---|
| 製品名 | Anti-Galectin 3 antibody |
| 製品の詳細 | Rabbit polyclonal to Galectin 3 |
| 特異性 | This antibody reacts specifically with 32 kDa human galectin-3 protein from chondrocytes and HT29 cell lines. |
| アプリケーション | 適用あり: IHC-P, ELISA, WB |
| 種交差性 | 交差種: Human |
| 免疫原 | 20 residue N-terminal synthetic peptide (Human). |

製品の特性

| | |
|--------|--|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
| バッファー | Preservative: None Constituents: Whole Serum |
| 精製度 | Whole antiserum |
| ポリ/モノ | ポリクローナル |
| アイソタイプ | IgG |

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab2473** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|--|
| IHC-P | | Use at an assay dependent concentration. |
| ELISA | | Use at an assay dependent concentration. |
| WB | | 1/200 - 1/5000. |

ターゲット情報

| | |
|-------|---|
| 機能 | Galactose-specific lectin which binds IgE. May mediate with the alpha-3, beta-1 integrin the stimulation by CSPG4 of endothelial cells migration. Together with DMBT1, required for terminal differentiation of columnar epithelial cells during early embryogenesis. |
| 組織特異性 | A major expression is found in the colonic epithelium. It is also abundant in the activated macrophages. |
| 配列類似性 | Contains 1 galectin domain. |
| 細胞内局在 | Nucleus. Cytoplasmic in adenomas and carcinomas. May be secreted by a non-classical secretory pathway and associate with the cell surface. |

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors