

Product datasheet

Glycophorin C overexpression 293T lysate (whole cell) ab94093

画像数 2

製品の概要

製品名	Glycophorin C overexpression 293T lysate (whole cell)
特記事項	ab94093 is a 293T cell transfected lysate in which Human Glycophorin C has been transiently over-expressed using a pCMV-Glycophorin C plasmid. The lysate is provided in 1X Sample Buffer. <b>Note:</b> For more details about how the transfected lysate was prepared <a href="#">view preparation notes</a>
アプリケーション	適用あり: WB

法規制情報

医薬用外劇物

製品の特性

製品の状態	Liquid
保存方法	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 6.80 Constituent: 100% 1x Sample Buffer
精製度	Whole Cell Lysate
背景	Function: This protein is a minor sialoglycoprotein in human erythrocyte membranes. The blood group Gerbich antigens and receptors for Plasmodium falciparum merozoites are most likely located within the extracellular domain. Glycophorin-C plays an important role in regulating the stability of red cells. Tissue specificity: Glycophorin-C is expressed in erythrocytes. Glycophorin-D is ubiquitous. Similarity: Belongs to the glycophorin-C family.


アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab94093** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent dilution.

## 画像



ab94093 at 15µg/lane on an SDS-PAGE gel.



**All lanes :** Anti-Glycophorin C antibody  
([ab55222](#)) at 1/500 dilution

**Lane 1 :** Glycophorin C overexpression 293T  
lysate (whole cell) (ab94093)

**Lane 2 :** 293T non-transfected lysate

Lysates/proteins at 25 µg per lane.

### Secondary

**All lanes :** Goat Anti-mouse IgG (H and L)  
HRP conjugated at 1/2500 dilution

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors