

## Product datasheet

# Goat Anti-Armenian hamster IgG H&L (Alkaline Phosphatase) ab5746

### 製品の概要

製品名	Goat Anti-Armenian hamster IgG H&L (Alkaline Phosphatase)
由来種	Goat
ターゲット生物種	Armenian hamster
特異性	Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Alkaline Phosphatase (calf intestine), anti-Goat Serum, Armenian Hamster and Armenian Hamster Serum.
アプリケーション	<b>適用あり:</b> ICC/IF, Dot blot, WB, ELISA, IHC
免疫原	Armenian Hamster IgG whole molecule
標識	Alkaline Phosphatase

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.1% Sodium Azide Constituents: 50% Glycerol, 10mg/ml BSA, 0.15M Sodium chloride, 0.05M Tris HCl, 1mM Magnesium chloride, 0.1mM Zinc chloride, pH 8
特記事項 (精製)	This product was prepared from monospecific antiserum by immunoaffinity chromatography using Armenian Hamster IgG coupled to agarose.
特記事項 (標識)	Alkaline Phosphatase (Calf Intestine) (Molecular Weight 140,000 daltons).
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab5746** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ICC/IF		1/200 - 1/1000.
Dot blot		Use at an assay dependent dilution.
WB		1/500 - 1/2500.
ELISA		1/2000 - 1/25000.
IHC		Use at an assay dependent dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors