

## Product datasheet

# Goat Anti-Cat IgG Fc (FITC) ab112800

### 製品の概要

製品名	Goat Anti-Cat IgG Fc (FITC)
製品の詳細	Goat polyclonal Secondary Antibody to Cat IgG - Fc (FITC)
由来種	Goat
ターゲット生物種	Cat
アプリケーション	<b>適用あり:</b> ICC/IF, IHC-P
標識	FITC. Ex: 493nm, Em: 528nm

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.09% Sodium azide Constituents: 99% PBS, 0.2% BSA
精製度	Immunogen affinity purified
特記事項(精製)	Antiserum was solid phase adsorbed to ensure class specificity. ab112800 was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to fluorescein isothiocyanate (FITC).
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab112800** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ICC/IF		1/200 - 1/2000.
IHC-P		1/200 - 1/2000.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors