

## Product datasheet

# Goat Anti-Chicken IgY H&L (Texas Red ®) ab6875

### 画像数 1

#### 製品の概要

製品名	Goat Anti-Chicken IgY H&L (Texas Red ®)
由来種	Goat
ターゲット生物種	Chicken
アプリケーション	<b>適用あり:</b> Immunomicroscopy, Flow Cyt, ICC/IF, IHC-P, IHC-Fr, ELISA
免疫原	chicken IgY whole molecule
標識	Texas Red ®. Ex: 596nm, Em: 620nm

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.60 Preservative: 0.01% Sodium azide Constituents: 0.88% Sodium chloride, 1% BSA, 0.42% Potassium phosphate
精製度	Affinity purified
特記事項 (精製)	This product was prepared from monospecific antiserum by immunoaffinity chromatography using chicken IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
特記事項 (標識)	Texas Red <sup>TM</sup> Sulfonyl Chloride (Molecular Weight 625 daltons) Absorption Wavelength: 596 nm Emission Wavelength: 620 nm Fluorochrome/Protein Ratio: 4.9 moles Texas Red <sup>ä</sup> per mole of Rabbit IgG
ポリ/モノ	ポリクローナル

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab6875** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Immunomicroscopy		

Flow Cyt

ICC/IF

IHC-P


IHC-Fr

ELISA

#### 追加情報

suitable for immunomicroscopy and flow cytometry or FACS analysis as well as other antibody based fluorescent assays requiring lot-to-lot consistency

#### 画像



Ab6875 was used at dilution 1/400 with the primary antibody [ab24525](#) in ICC/IF. See the review on [ab24525](#).

Immunocytochemistry/ Immunofluorescence - Goat polyclonal Secondary Antibody to Chicken IgY - H&L (TR) (ab6875)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors