

## Product datasheet

# Goat Anti-Guinea pig IgG H&L ab97151

### 製品の概要

製品名	Goat Anti-Guinea pig IgG H&L
ターゲット生物種	Guinea pig
特異性	By immunoelectrophoresis and ELISA this antibody reacts specifically with Guinea Pig IgG and with light chains common to other Guinea Pig immunoglobulins. No antibody was detected against non immunoglobulin serum proteins.
アプリケーション	<b>適用あり:</b> ICC, IHC-P, ELISA, WB

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.09% Sodium Azide Constituents: PBS
精製度	Immunogen affinity purified
特記事項 (精製)	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab97151** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ICC		Use at an assay dependent dilution.
IHC-P		1/200 - 1/2000.
ELISA		1/100 - 1/30000. ELISA-Coating 1/100 - 1/500; ELISA-Primary 1/1000 - 1/30000
WB		1/1000 - 1/30000. WB-Colorimetric 1/1000 - 1/10000; WB-Chemiluminescent 1/1000 - 1/30000

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors