

Product datasheet

Goat Anti-Human IgA alpha chain (Alkaline Phosphatase) ab97212

製品の概要

| | |
|----------|--|
| 製品名 | Goat Anti-Human IgA alpha chain (Alkaline Phosphatase) |
| ターゲット生物種 | Human |
| 特異性 | By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgA. Cross reactivity with other immunoglobulins and light chains is less than 0.1%. |
| アプリケーション | 適用あり: ICC, IHC-P, ELISA, WB |
| 標識 | Alkaline Phosphatase |

製品の特性

| | |
|----------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C. |
| バッファー | Preservative: 0.09% Sodium Azide Constituents: 0.2% BSA, 0.1M Sodium chloride, 50mM HEPES, 1mM Magnesium chloride, 0.1mM Zinc chloride, pH 7.1 |
| 精製度 | Immunogen affinity purified |
| 特記事項(精製) | Antiserum was solid phase adsorbed to ensure class specificity. This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to Alkaline Phosphatase (AP). |
| ポリ/モノ | ポリクローナル |
| アイソタイプ | IgG |

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab97212** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|-------------------------------------|
| ICC | | Use at an assay dependent dilution. |
| IHC-P | | 1/100 - 1/1000. |

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|--------------------------------------|
| ELISA | | 1/1000 - 1/10000. (Primary ELISA) |
| WB | | 1/1000 - 1/10000. (Chemiluminescent) |

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors