

## Product datasheet

# Goat IgG (Fc) - Isotype Control ab37388

### 3 References

#### 製品の概要

|          |  |
|----------|--|
| 製品名      | Goat IgG (Fc) - Isotype Control  |
| 特異性      | This antibody has been selected to be an isotype control as it has no known specificity. |
| アプリケーション | <b>適用あり:</b> ELISA   |

#### 製品の特性

|        |  |
|--------|--|
| 製品の状態  | Liquid   |
| 保存方法   | Shipped at 4°C. Store at +4°C.                       |
| バッファー  | pH: 8.20<br>Constituent: 100% Borate buffered saline |
| 精製度    | IgG fraction   |
| ポリ/モノ  | ポリクローナル  |
| アイソタイプ | IgG  |

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab37388** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項                                |
|----------|-----------|-------------------------------------|
| ELISA    |           | Use at an assay dependent dilution. |

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors