

## Product datasheet

# Goat Anti-Mouse IgG F(ab) (TRITC) preadsorbed ab50596

1 References [画像数 1](#)

### 製品の概要

|              |   |                              |
|--------------|---|------------------------------|
| 製品名          | Goat Anti-Mouse IgG F(ab) (TRITC) preadsorbed   |                              |
| ターゲット生物種     | Mouse   |                              |
| 特異性          | The conjugate can detect mouse IgA, mouse IgM, and all subclasses of mouse IgG. The conjugate shows no reactivity with mouse IgG Fc fragment. |                              |
| アプリケーション     | <b>適用あり:</b> IHC-P, ICC/IF  |                              |
| 吸着処理血清       | Human   | <a href="#">more details</a> |
| 免疫原          | Purified mouse IgG, Fab fragment  |                              |
| ポジティブ・コントロール | Human peripheral blood lymphocytes and human tonsil sections  |                              |
| 標識           | TRITC. Ex: 547nm, Em: 572nm   |                              |

### 製品の特性

|           |  |
|-----------|--|
| 製品の状態     | Liquid   |
| 保存方法      | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.                         |
| バッファー     | Preservative: 15mM Sodium Azide<br>Constituents: 0.01M PBS, pH 7.4   |
| 精製度       | Immunogen affinity purified  |
| 特記事項 (精製) | The antibody preparation is solid phase adsorbed with human IgG to ensure minimal cross reactivity in tissue or cell preparations.                       |
| 特記事項 (標識) | Ab50596 is conjugated to crystalline tetramethylrhodamine isothiocyanate (TRITC) in an alkaline reaction, then further purified to remove unbound TRITC. |
| ポリ/モノ     | ポリクローナル  |
| アイソタイプ    | IgG  |


### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab50596** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項   |
|----------|-----------|--------|
| IHC-P    |           | 1/100. |
| ICC/IF   |           | 1/32.  |

## 画像



Immunocytochemistry/ Immunofluorescence - IgG antibody (TRITC) (ab50596)

Immunofluorescent staining of microglia in the hippocampus (CA3 region) in rat brain. Rat brain sections stained with Monoclonal Anti-Glial Fibrillary Acidic Protein (clone no. G-A-5), and ab50596 at 1/32 dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors