

Product datasheet

Goat Anti-Mouse IgG H&L preadsorbed ab98800

製品の概要

| | |
|----------|---|
| 製品名 | Goat Anti-Mouse IgG H&L preadsorbed |
| 製品の詳細 | Goat polyclonal Secondary Antibody to Mouse IgG - H&L, pre-adsorbed |
| ターゲット生物種 | Mouse |
| アプリケーション | 適用あり: IHC-P, ELISA, WB, ICC |
| 吸着処理血清 | Human, Rat more details |

製品の特性

| | |
|-----------|--|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
| バッファー | Preservative: 0.09% Sodium Azide Constituents: PBS |
| 精製度 | Immunogen affinity purified |
| 特記事項 (精製) | This antibody was isolated by affinity chromatography using antigen coupled to agarose beads. |
| ポリモノ | ポリクローナル |
| アイソタイプ | IgG |

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab98800** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|---|
| IHC-P | | 1/200 - 1/2000. |
| ELISA | | 1/100 - 1/30000. Coating: Use 1/100-1/500 dilution; Primary: Use 1/1,000-1/30,000 dilution. |
| WB | | 1/1000 - 1/30000. Colorimetric: Use 1/1,000-1/10,000 dilution; Chemiluminescent: Use 1/1,000-1/30,000 dilution. |

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|-------------------------------------|
| ICC | | Use at an assay dependent dilution. |

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors