

Product datasheet

Goat Anti-Mouse IgG+IgM H&L (12nm Gold) preadsorbed ab105289

1 References

製品の概要

製品名	Goat Anti-Mouse IgG+IgM H&L (12nm Gold) preadsorbed	
由来種	Goat	
ターゲット生物種	Mouse	
アプリケーション	適用あり: ELISA, Electron Microscopy, WB	
吸着処理血清	Cow, Horse, Human	more details
標識	Gold 12nm	

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 8.50 Preservative: 0.05% Sodium azide Constituents: 0.381% Sodium borate, 1.5% BSA, 0.164% Sodium phosphate, 0.87% Sodium chloride
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab105289** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent concentration.
Electron Microscopy		1/20 - 1/40.

アプリケーション

Abreviews

特記事項

WB

Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors