

Product datasheet

Goat Anti-Mouse IgG+IgM H&L (Texas Red ®) preadsorbed ab47845

製品の概要

製品名	Goat Anti-Mouse IgG+IgM H&L (Texas Red ®) preadsorbed
由来種	Goat
ターゲット生物種	Mouse
特異性	This antibody reacts with Mouse heavy and light chains (IgG and IgM). Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Goat Serum, Mouse IgG and Mouse IgM.
アプリケーション	適用あり: Flow Cyt, ICC/IF, ELISA
吸着処理血清	Cow, Horse, Human more details
免疫原	Full length Mouse IgG and IgM whole molecules.
標識	Texas Red ®. Ex: 596nm, Em: 620nm

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.01% Sodium Azide Constituents: 10mg/ml BSA, 0.15M Sodium chloride, 0.02M Potassium phosphate, pH 7.2
精製度	Immunogen affinity purified
特記事項 (精製)	Prepared from polyspecific antiserum by immunoaffinity chromatography using antigens coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities.
特記事項 (標識)	Label: Texas Red Sulfonyl Chloride (Molecular Weight 625 daltons) Absorption Wavelength: 596 nm : Emission Wavelength: 620 nm 3.7 moles Texas Red per mole of Goat IgG
ポリモノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab47845** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Flow Cyt		1/500 - 1/2500.
ICC/IF		1/1000 - 1/5000.
ELISA		1/10000 - 1/50000.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors