

Product datasheet

Goat Anti-Mouse IgM (SureLight® Allophycocyanin) ab99905

製品の概要

製品名	Goat Anti-Mouse IgM (SureLight® Allophycocyanin)
由来種	Goat
ターゲット生物種	Mouse
アプリケーション	適用あり: Flow Cyt
標識	SureLight® Allophycocyanin. Ex: 652nm, Em: 657nm

製品の特性

製品の状態	Lyophilised: Lyophilized material is stable for up to one year. After product has been reconstituted, product should be stored at 2-8°C in the dark and be used within 3 months. If further dilution of the conjugate is required, use diluted material within one week.
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: None Constituents: 0.1% BSA, 50mM Sucrose, 150mM Sodium chloride, 1 ug/ml pentachlorophenol pH 8.2
精製度	Protein G purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab99905** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Flow Cyt		Use a concentration of 3 µg/ml.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors