

Product datasheet

Goat Anti-Mouse IgM mu chain (6nm Gold) ab39606

製品の概要

製品名	Goat Anti-Mouse IgM mu chain (6nm Gold)
由来種	Goat
ターゲット生物種	Mouse
アプリケーション	適用あり: Electron Microscopy
免疫原	Full length protein (Mouse)
標識	Gold 6nm

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 15mM Sodium Azide Constituents: 1% BSA, PBS
精製度	Immunogen affinity purified
特記事項 (精製)	This product has been immuno-affinity purified and immuno cross-absorbed to reduce non-specific reactions. The activity of each lot is determined using a dot-spot test system as described by Moeremans et al., J. Immunol. Methods, 74, (1984), 353.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG
特記事項	This product is an immunogold reagent built around colloidal gold particles of 6nm. The particle population is monodisperse and thus shows minimal size variation and overlap. Typically, the coefficient of variance for the 6nm particle size conjugate is less than 12 %.

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab39606** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Electron Microscopy		

追加情報

Electron Microscopy: Use at an assay dependent dilution. The method of production of this product allows the small particle size of 6nm instead of 5nm with the same sensitivity but increased visibility in the standard transmission electron microscope (see also van Bergen et en Henegouwen & Leunissen, Histochemistry, 85, (1986), 81.)

Light microscopy: Use at an assay dependent dilution.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors