

Product datasheet

Goat Anti-Rabbit IgG Fc (DyLight® 550) preadsorbed ab98463

画像数 1

製品の概要

製品名	Goat Anti-Rabbit IgG Fc (DyLight® 550) preadsorbed	
ターゲット生物種	Rabbit	
アプリケーション	適用あり: Flow Cyt, ICC/IF, IHC-P	
吸着処理血清	Human, Mouse, Rat	more details
標識	DyLight® 550. Ex: 562nm, Em: 576nm	

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.09% Sodium Azide Constituents: 0.2% BSA, PBS
精製度	Immunogen affinity purified
特記事項 (精製)	This antibody was isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to DyLight® 550.
ポリモノ	ポリクローナル
アイソタイプ	IgG
特記事項	DyLight® 550 replaces DyLight® 549.


アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab98463** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Flow Cyt		1/50 - 1/200.
ICC/IF		1/50 - 1/500.
IHC-P		1/50 - 1/500.

画像



DyLight®-Goat polyclonal Secondary Antibody to Rabbit IgG - Fc (DyLight® 550), pre-adsorbed(ab98463)

Emission spectra of DyLight® fluorochromes available in our catalog.

Line colors represent the approximate visible colors of the wavelength maxima.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors