

## Product datasheet

# Goat Anti-Rabbit IgG H&L (10nm Gold) preadsorbed ab27234

### 9 References

#### 製品の概要

製品名	Goat Anti-Rabbit IgG H&L (10nm Gold) preadsorbed	
由来種	Goat	
ターゲット生物種	Rabbit	
アプリケーション	<b>適用あり:</b> Electron Microscopy	
吸着処理血清	Human	<a href="#">more details</a>
免疫原	Full lenght rabbit IgG (H and L).	
標識	Gold 10nm	

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 15mM Sodium Azide Constituents: 20% Glycerol, 1% BSA, Tris buffered saline. pH 8.2
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab27234** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Electron Microscopy		Use at an assay dependent dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors