

## Product datasheet

# Goat Anti-Turkey IgY H&L (Alkaline Phosphatase) ab125855

### 製品の概要

製品名	Goat Anti-Turkey IgY H&L (Alkaline Phosphatase)
製品の詳細	Goat polyclonal Secondary Antibody to Turkey IgY - H&L (AP)
由来種	Goat
ターゲット生物種	Turkey
アプリケーション	<b>適用あり:</b> ICC/IF, ELISA, WB, Dot blot, IHC-P, Other
免疫原	The details of the immunogen for this antibody are not available.
標識	Alkaline Phosphatase

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 8.00 Preservative: 0.1% Sodium azide Constituents: 0.01% Magnesium chloride, 0.61% Tris, 50% Glycerol
精製度	IgG fraction
特記事項 (精製)	Affinity chromatography on Chicken IgG covalently linked to agarose.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG
特記事項	Pooled antisera from goats hyperimmunized with Turkey and Chicken IgG.

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab125855** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ICC/IF		Use at an assay dependent concentration.
ELISA		1/2000 - 1/4000.

アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent concentration.
Dot blot		Use at an assay dependent concentration.
IHC-P		Use at an assay dependent concentration.
Other		Use at an assay dependent concentration. Fluorescent-Linked-Immunosorbent-Assay (FLISA)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors