

Product datasheet

Anti-HSC70 Interacting Protein HIP antibody ab117109

画像数 1

製品の概要

製品名	Anti-HSC70 Interacting Protein HIP antibody
製品の詳細	Rabbit polyclonal to HSC70 Interacting Protein HIP
由来種	Rabbit
アプリケーション	<b>適用あり:</b> IHC-P, ELISA
種交差性	<b>交差種:</b> Human
免疫原	Synthetic peptide: PDSKKVEEDLKADEPS conjugated to KLH, corresponding to amino acids 60-75 of Human HSC70 Interacting Protein HIP (NP_003923.2). <a href="#">Run BLAST with</a> <a href="#">Run BLAST with</a>
ポジティブ・コントロール	Human adrenal tissue.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at -20°C.
バッファー	Preservative: 0.1% Sodium azide Constituent: 99% PBS
精製度	Protein G purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab117109** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項

アプリケーション	Abreviews	特記事項
IHC-P		Use a concentration of 10 µg/ml. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
ELISA		1/1 - 1/1000.

## ターゲット情報

機能	One HIP oligomer binds the ATPase domains of at least two HSC70 molecules dependent on activation of the HSC70 ATPase by HSP40. Stabilizes the ADP state of HSC70 that has a high affinity for substrate protein. Through its own chaperone activity, it may contribute to the interaction of HSC70 with various target proteins.
配列類似性	Belongs to the FAM10 family. Contains 1 ST11 domain. Contains 3 TPR repeats.
細胞内局在	Cytoplasm.

## 画像



ab117109 at 10µg/ml staining HSC70 Interacting Protein HIP in Formalin-fixed, Paraffin-embedded Human Adrenal tissue by Immunohistochemistry followed by biotinylated secondary antibody, alkaline phosphatase-streptavidin and chromogen.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-HSC70 Interacting Protein HIP antibody (ab117109)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors