

Product datasheet

Anti-HSPA4L antibody ab81221

★★★★☆ 3 Abreviews 1 References 画像数 1

製品の概要

製品名	Anti-HSPA4L antibody
製品の詳細	Rabbit polyclonal to HSPA4L
由来種	Rabbit
アプリケーション	適用あり: ELISA, WB
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Rabbit, Horse, Chicken, Guinea pig, Cow, Cat, Dog, Pig 
免疫原	Synthetic peptide, corresponding to a region within C terminal amino acids 720-769 (RNKDERYDHLDPTEMEKVEKICISDAMSWLNSKMNAQNKLSLTQDPVVKVS) of human HSPA4L (NP_055093). Run BLAST with Run BLAST with
ポジティブ・コントロール	Fetal brain lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.09% Sodium azide Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
特記事項(精製)	Purified by peptide affinity chromatography method.
ポリモノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab81221** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent concentration. ELISA titre using peptide based assay: 1:312500
WB	★★★★☆	Use a concentration of 1 µg/ml. Predicted molecular weight: 94 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

ターゲット情報

関連性	The HSPA4L gene, also known as Apg1 or Osp94, belongs to the HSP110 heat shock gene family, which includes three genes encoding highly conserved proteins. HSPA4L is expressed ubiquitously and predominantly in the testis.
細胞内局在	Cytoplasm; Nucleus. Note= May translocate to the nucleus after heat shock

画像



Anti-HSPA4L antibody (ab81221) at 1 µg/ml
(in 5% skim milk / PBS buffer) + fetal brain
lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000
dilution

Predicted band size: 94 kDa

Observed band size: 90 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors