

Product datasheet

Anti-HSV1 antibody, prediluted ab933

製品の概要

製品名	Anti-HSV1 antibody, prediluted
製品の詳細	Rabbit polyclonal to HSV1, prediluted
アプリケーション	適用あり: IHC-P, IHC-Fr
免疫原	Whole rabbit corneal cells, infected with herpes simplex virus type 1 (strain Mac Intyre), were solubilized with detergent.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C. Do Not Freeze.
バッファー	Preservative: 0.05% Sodium Azide Constituents: Inert stabilizer
精製度	Protein A purified
ポリ/モノ	ポリクローナル
ミエローマ	unknown
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab933** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
IHC-P		
IHC-Fr		

追加情報

This antibody has been pretitered and quality controlled to work on formalin-fixed paraffin-embedded as well as acetone fixed cryostat tissue sections.
No further titration is required.
We suggest an incubation period of 30-60 minutes at room temperature.
However, depending upon the fixation conditions and the staining system employed, optimal

incubation conditions should be determined by the user.

ターゲット情報

関連性 Herpes simplex type 1 (HSV-1) belongs to a family that includes HSV-2, Epstein-Barr virus (EBV) and Varicella zoster (chicken pox) virus amongst others. HSV-1 and HSV-2 are extremely difficult to distinguish from each other. Members of this family have a characteristic virion structure. The double stranded DNA genome is contained within an icosahedral capsid embedded in a proteinaceous layer (tegument) and surrounded by a lipid envelope, derived from the nuclear membrane of the last host, which is decorated with virus-specific glycoproteins spikes. These viruses are capable of entering a latent phase where the host shows no visible sign of infection and levels of infectious agent become very low. During the latent phase the viral DNA is integrated into the genome of the host cell.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors