

Product datasheet

Human APXL peptide ab23127

製品の概要

製品名 Human APXL peptide

製品の詳細

由来 Synthetic

アミノ酸配列

生物種 Human

配列 KCLLDSLQPERGK

領域 1604 to 1616

特性

Our [Abpromise guarantee](#) covers the use of **ab23127** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション Blocking

製品の状態 Liquid

前処理および保存

保存方法および安定性 Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

関連情報

機能 May be involved in endothelial cell morphology changes during cell spreading. In the retinal pigment epithelium, may regulate the biogenesis of melanosomes and promote their association with the apical cell surface by inducing gamma-tubulin redistribution.

組織特異性 Abundant in retina and melanoma; also in brain, placenta, lung, kidney and pancreas.

配列類似性 Belongs to the shroom family.
Contains 1 ASD1 domain.
Contains 1 ASD2 domain.
Contains 1 PDZ (DHR) domain.

ドメイン	The ASD1 domain mediates F-actin binding.
細胞内局在	Apical cell membrane. Cell junction > tight junction. Cytoplasm > cytoskeleton. Associates with cortical F-actin.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors