

Product datasheet

Human DAGLB peptide ab104607

製品の概要

製品名 Human DAGLB peptide

製品の詳細

由来 Synthetic

アミノ酸配列

生物種 Human

配列 STAELFSTYF SDTDLVPSDI AAGLALLHQQ
QDNIRNNQEP AQQVCHAPGS

領域 216 to 265

特性

Our [Abpromise guarantee](#) covers the use of **ab104607** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション Blocking - Blocking peptide for Anti-DAGLB antibody ([ab84733](#))

製品の状態 Liquid

前処理および保存

保存方法および安定性 Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

Preservative: None

Constituents: dH2O

関連情報

機能 Catalyzes the hydrolysis of diacylglycerol (DAG) to 2-arachidonoyl-glycerol (2-AG), the most abundant endocannabinoid in tissues. Required for axonal growth during development and for retrograde synaptic signaling at mature synapses.

配列類似性 Belongs to the AB hydrolase superfamily. Lipase family.

細胞内局在 Cell membrane.

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors