

Product datasheet

Human DFFB peptide ab8439

製品の概要

製品名 Human DFFB peptide

製品の詳細

由来 Synthetic

アミノ酸配列

生物種 Human

配列 NGSYFDRGAKGGSRLC Amino acids 203 to 218 of human DFF40 (1,2)

領域 203 to 218

特性

Our [Abpromise guarantee](#) covers the use of **ab8439** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション Blocking
Western blot

製品の状態 Liquid

前処理および保存

保存方法および安定性 Shipped at 4°C. After reconstitution store at -20°C. Avoid freeze / thaw cycles.
PBS with 0.1% BSA 0.02% sodium azide pH7.2

関連情報

機能 Nuclease that induces DNA fragmentation and chromatin condensation during apoptosis.
Degrades naked DNA and induces apoptotic morphology.

配列類似性 Contains 1 CIDE-N domain.

細胞内局在 Cytoplasm. Nucleus.

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors