

## Product datasheet

# Human DOK1 peptide ab39756

### 製品の概要

製品名 Human DOK1 peptide

### 製品の詳細

由来 Synthetic

#### アミノ酸配列

生物種 Human

配列 PQGPAFPEPGTATGS

領域 425 to 439

### 特性

Our [Abpromise guarantee](#) covers the use of **ab39756** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション Western blot  
Blocking - Blocking peptide for Anti-DOK1 antibody ([ab8112](#))

製品の状態 Liquid

### 前処理および保存

保存方法および安定性 Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.  
Preservative: 0.02% Sodium Azide  
Constituents: 0.01% BSA, PBS, pH 7.2

### 関連情報

機能 DOK proteins are enzymatically inert adaptor or scaffolding proteins. They provide a docking platform for the assembly of multimolecular signaling complexes. DOK1 appears to be a negative regulator of the insulin signaling pathway. Modulates integrin activation by competing with talin for the same binding site on ITGB3.

組織特異性 Expressed in pancreas, heart, leukocyte and spleen. Expressed in both resting and activated peripheral blood T-cells.

<b>配列類似性</b>	Belongs to the DOK family. Type A subfamily. Contains 1 IRS-type PTB domain. Contains 1 PH domain.
<b>ドメイン</b>	The PTB domain mediates receptor interaction.
<b>翻訳後修飾</b>	Constitutively tyrosine-phosphorylated. Phosphorylated on tyrosine residues by the insulin receptor kinase. Results in the negative regulation of the insulin signaling pathway.
<b>細胞内局在</b>	Cytoplasm and Cytoplasm > perinuclear region.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors