

Product datasheet

Human FOXN1 peptide ab45510

製品の概要

製品名 Human FOXN1 peptide

製品の詳細

由来 Synthetic

アミノ酸配列

生物種 Human

配列 C-SVYLSPSSKPVALA

領域 635 to 648

特性

Our [Abpromise guarantee](#) covers the use of **ab45510** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション Blocking - Blocking peptide for Anti-FOXN1 antibody ([ab26069](#))

製品の状態 Liquid

前処理および保存

保存方法および安定性 Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

関連情報

機能 Transcriptional regulator involved in development.

組織特異性 Expressed in thymus.

関連疾患 Defects in FOXN1 are the cause of T-cell immunodeficiency congenital alopecia and nail dystrophy (TIDAND) [MIM:601705]. A disorder characterized by the association of congenital alopecia, severe T-cell immunodeficiency, and ridging and pitting of all nails.

配列類似性 Contains 1 fork-head DNA-binding domain.

細胞内局在 Nucleus.

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors