

## Product datasheet

# Human MSH4 peptide ab23238

### 製品の概要

製品名 Human MSH4 peptide

### 製品の詳細

由来 Synthetic

#### アミノ酸配列

生物種 Human

配列 C-FPRTEQVPEKTEE

領域 924 to 936

### 特性

Our [Abpromise guarantee](#) covers the use of **ab23238** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション Blocking

製品の状態 Liquid

### 前処理および保存

保存方法および安定性 Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

### 関連情報

機能 Involved in meiotic recombination. Required for reciprocal recombination and proper segregation of homologous chromosomes at meiosis.

組織特異性 Testis and ovary specific.

配列類似性 Belongs to the DNA mismatch repair MutS family.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors