

Product datasheet

Human PITPN peptide ab23201

製品の概要

製品名 Human PITPN peptide

製品の詳細

由来 Synthetic

アミノ酸配列

生物種 Human

配列 C-RQKDPVKGMTADD

領域 258 to 270

特性

Our [Abpromise guarantee](#) covers the use of **ab23201** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション Blocking - Blocking peptide for Anti-PITPN antibody ([ab10265](#))

製品の状態 Liquid

前処理および保存

保存方法および安定性 Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

関連情報

機能 Catalyzes the transfer of PtdIns and phosphatidylcholine between membranes.

組織特異性 Expressed in a wide range of tissues.

配列類似性 Belongs to the PtdIns transfer protein family. PI transfer class I subfamily.

細胞内局在 Cytoplasm.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors