

Product datasheet

Human PLRG1 peptide ab22860

製品の概要

製品名 Human PLRG1 peptide

製品の詳細

由来 Synthetic

アミノ酸配列

生物種 Human

配列 C-PVSWKPEIIKRRKF

領域 501 to 514

特性

Our [Abpromise guarantee](#) covers the use of **ab22860** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション Blocking - Blocking peptide for Anti-PLRG1 antibody ([ab1372](#))

製品の状態 Liquid

前処理および保存

保存方法および安定性 Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

関連情報

機能 Necessary for spliceosome assembly and for pre-mRNA splicing.

配列類似性 Belongs to the WD repeat PRL1/PRL2 family.
Contains 7 WD repeats.

細胞内局在 Nucleus speckle.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors