

Product datasheet

Human PSCDBP peptide ab22893

製品の概要

製品名 Human PSCDBP peptide

製品の詳細

由来 Synthetic

アミノ酸配列

生物種 Human

配列 SLQRLLQHSSNGN-C

領域 2 to 14

特性

Our [Abpromise guarantee](#) covers the use of **ab22893** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション Blocking - Blocking peptide for Anti-PSCDBP antibody ([ab2247](#))

製品の状態 Liquid

前処理および保存

保存方法および安定性 Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

関連情報

機能 By its binding to cytohesin-1 (CYTH1), it modifies activation of ARFs by CYTH1 and its precise function may be to sequester CYTH1 in the cytoplasm.

組織特異性 Expressed in lymph nodes, thymus, spleen, lung, peripheral blood leukocytes and bone marrow.

配列類似性 Contains 1 PDZ (DHR) domain.

細胞内局在 Cytoplasm. Early endosome. Recruited from the cytosol to endosomes by SNX27.

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors