

## Product datasheet

# Human TGN46 peptide ab41822

### 製品の概要

製品名 Human TGN46 peptide

### 製品の詳細

由来 Synthetic

#### アミノ酸配列

生物種 Human

配列 VVPEQPSWKDHSKPIS

領域 252 to 267

### 特性

Our [Abpromise guarantee](#) covers the use of **ab41822** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション Neutralising  
Blocking - Blocking peptide for Anti-TGN46 antibody ([ab41702](#))

製品の状態 Liquid

備考 Purity of peptide: >70%

### 前処理および保存

保存方法および安定性 Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

### 関連情報

機能 May be involved in regulating membrane traffic to and from trans-Golgi network.

組織特異性 Isoform TGN46 is widely expressed. Isoform TGN51 is more abundant in fetal lung and kidney. Isoform TGN48 is barely expressed in embryonic kidney and promyelocytic cells.

細胞内局在 Cell membrane. Golgi apparatus > trans-Golgi network membrane. Primarily in trans-Golgi network. Cycles between the trans-Golgi network and the cell surface returning via endosomes.

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors