

Product datasheet

Human Wnt3 peptide ab50436

製品の概要

製品名 Human Wnt3 peptide

製品の詳細

由来 Synthetic

アミノ酸配列

生物種 Human

配列 CGRGHNTRTEKRKEK

領域 315 to 329

特性

Our [Abpromise guarantee](#) covers the use of **ab50436** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション Blocking - Blocking peptide for Anti-Wnt3 antibody ([ab50341](#))

製品の状態 Liquid

前処理および保存

保存方法および安定性 Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

関連情報

機能 Ligand for members of the frizzled family of seven transmembrane receptors. Wnt-3 and Wnt-3a play distinct roles in cell-cell signaling during morphogenesis of the developing neural tube.

関連疾患 Defects in WNT3 are the cause of autosomal recessive tetra-amelia (ARTTRA) [MM:273395]. Tetra-amelia is a rare human genetic disorder characterized by complete absence of all four limbs and other anomalies such as craniofacial, nervous system, pulmonary, skeletal and urogenital defects.

配列類似性 Belongs to the Wnt family.

細胞内局在 Secreted > extracellular space > extracellular matrix.

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors