

## Product datasheet

### Human ZNF179 peptide ab23057

#### 製品の概要

---

製品名 Human ZNF179 peptide

#### 製品の詳細

---

由来 Synthetic

#### アミノ酸配列

生物種 Human

配列 C-GDREPLLQEE

領域 622 to 632

#### 特性

---

Our [Abpromise guarantee](#) covers the use of **ab23057** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション Blocking - Blocking peptide for Anti-ZNF179 antibody ([ab4838](#))

製品の状態 Liquid

#### 前処理および保存

---

保存方法および安定性 Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

#### 関連情報

---

関連性 ZNF179 encodes a member of the RING finger protein family of transcription factors. The protein is primarily expressed in brain. The gene is located within the Smith-Magenis syndrome region on chromosome 17.

細胞内局在 Cytoplasmic, Nuclear and integral to membrane.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors